


Rabbit Anti-Human Cd226/Dnam-1 Monoclonal Antibody For Flow Cytometry - Q15762

Item Number: CM0007986



Introduction

High-performance -conjugated rabbit monoclonal antibody targeting human CD226/DNAM-1, validated for flow cytometry. Essential tool for studying NK and T-cell-mediated cytotoxicity in immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD226/DNAM-1 mAb
Remarks/Alias	PTA1; DNAM1; DNAM-1; TLISA1
Species	Human
GeneID (Human)	10666
GeneID	10666
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-247 of human CD226/DNAM-1 (NP_006557.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15762
Protein Weight	39kDa
Shipping	Ice bag