

Rabbit Anti-Human Cd228/Mfi2 Monoclonal Antibody For Flow Cytometry - P08582

Item Number: CM0012084



Introduction

-conjugated rabbit monoclonal anti-human CD228/MFI2 (Melanotransferrin) antibody for flow cytometry. Detects a 32/80 kDa protein in human samples. Ideal for melanoma, fetal tissue, and iron transport research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD228/MFI2 mAb
Remarks/Alias	MTf; MFI2; MTF1; CD228; MAP97
Species	Human
GeneID (Human)	4241
GeneID	4241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-708 of human CD228/MFI2 (NP_005920.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08582
Protein Weight	32kDa/80kDa
Shipping	Ice bag