


Rabbit Anti-Human Cd235A/Glycophorin A Mab For Fc - P02724

Item Number: CM0008434



Introduction

-conjugated rabbit monoclonal antibody targeting human CD235a/Glycophorin A. Validated for flow cytometry, reacts with human. Major erythrocyte membrane protein, MN blood group antigen. Ideal for blood research and infectious disease studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD235a/Glycophorin A mAb
Remarks/Alias	MN; GPA; MNS; GPSAT; PAS-2; CD235a; GP _E rik; HGpMiV; HGpMiXI; HGpSta(C)
Species	Human
GeneID (Human)	2993
GeneID	2993
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-91 of human CD235a (NP_002090.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02724
Protein Weight	13kDa/16kDa
Shipping	Ice bag