


Rabbit Anti-Human Cd235A/Glycophorin A Monoclonal Antibody For Fc - P02724

Item Number: CM0008475



Introduction

-conjugated rabbit monoclonal antibody targeting human CD235a/Glycophorin A for flow cytometry. Validated in human samples. Major erythrocyte membrane protein with MN blood group antigens.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD235a/Glycophorin A mAb
Remarks/Alias	MN; GPA; MNS; GPSAT; PAS-2; CD235a; GPErik; HGpMiV; HGpMiXI; HGpSta(C)
Species	Human
Gene ID (Human)	2993
Gene ID	2993
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02724
Protein Weight	13kDa/16kDa
Shipping	Ice bag