


Rabbit Anti-Human Cd235A/Glycophorin A Monoclonal Antibody For Fc - P02724

Item Number: CM0011533



Introduction

-labeled rabbit monoclonal antibody targeting human CD235a/Glycophorin A (GYPA) for flow cytometry. Cross-reacts with human samples; suitable for erythrocyte membrane and blood group antigen research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD235a/Glycophorin A mAb
Remarks/Alias	MN; GPA; MNS; GPSAT; PAS-2; CD235a; GP _E rik; HGpMiV; HGpMiXI; HGpSta(C)
Species	Human
GeneID (Human)	2993
GeneID	2993
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 20-91 of human CD235a/Glycophorin A (NP_002090.4))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02724
Protein Weight	13kDa/16kDa
Shipping	Ice bag