


Rabbit Anti-Human Cd235A/Glycophorin A Monoclonal Antibody For Fc - P02724

Item Number: CM0005385



Introduction

High-quality rabbit monoclonal antibody specific for human Glycophorin A (CD235a), validated for flow cytometry. Ideal for erythrocyte membrane research, blood group antigen studies, and host-pathogen interaction assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD235a/Glycophorin A mAb
Remarks/Alias	MN, GPA, MNS, GPSAT, PAS-2, CD235a, GP _E rik, HGpMiV, HGpMiXI, HGpSta(C)
Species	Human
Gene ID	2993
Immunogen	Recombinant fusion protein containing amino acids 20-91 of human CD235a/Glycophorin A (NP_002090.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02724
Protein Weight	13kDa/16kDa
Shipping	Ice bag