

Rabbit Anti-Human Cd239/Bcam Monoclonal Antibody For Flow Cytometry - P50895

Item Number: CM0009370



Introduction

-conjugated rabbit monoclonal antibody targeting human CD239/BCAM. Validated for flow cytometry (FC), reacts specifically with human CD239, a 67kDa cell adhesion molecule involved in laminin binding and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD239/BCAM mAb
Remarks/Alias	BCAM; AU; CD239; LU; MSK19; basal cell adhesion molecule
Species	Human
GeneID (human)	4059
GeneID	4059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-547 of human CD239/BCAM (NP_005572.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P50895
Protein Weight	67kDa
Shipping	Ice bag