

## Rabbit Anti-Human Cd239/Bcam Monoclonal Antibody For Flow Cytometry - P50895

Item Number: CM0008091



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD239/BCAM, validated for flow cytometry. Ideal for studying cell adhesion, migration, and Lutheran blood group antigen in cancer and immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD239/BCAM mAb                                                                                       |
| Remarks/Alias    | BCAM                                                                                                                   |
| Species          | Human                                                                                                                  |
| GeneID (Human)   | 4059                                                                                                                   |
| GeneID           | 4059                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 32-547 of human CD239/BCAM (NP_005572.2) |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | FC                                                                                                                     |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | P50895                                                                                                                 |
| Protein Weight   | 67kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |