

Rabbit Anti-Human Cd24 Mab For If/Icc, Fc - P25063

Item Number: CM0005305

Introduction

Recombinant rabbit monoclonal antibody against human CD24, directly conjugated to . Validated for immunofluorescence and flow cytometry. Targets a GPI-anchored glycoprotein regulating B-cell activation and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD24 mAb
Remarks/Alias	CD24A; CD24
Species	Human
Gene ID	100133941
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P25063
Protein Weight	8kDa; 13kDa
Shipping	Ice bag