

Rabbit Anti-Human Cd24 Monoclonal Antibody For Fc - P25063

Item Number: CM0008637

Introduction

Rabbit anti-Human CD24 monoclonal antibody conjugated to for flow cytometry. Validated for FC, recognizing human CD24 involved in B-cell activation and cell differentiation. 8kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD24 mAb
Remarks/Alias	CD24A
Species	Human
GeneID (Human)	100133941
GeneID	100133941
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25063
Protein Weight	8kDa
Shipping	Ice bag