

Rabbit Anti-Human Cd244/2B4 Monoclonal Antibody For Flow Cytometry - Q9Bzw8

Item Number: CM0009571

Introduction

-conjugated rabbit monoclonal antibody targeting human CD244/2B4 (SLAMF4) for flow cytometry. Specific for natural killer and T cell subsets; validated in human samples. Ideal for immune cell profiling and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD244/2B4 mAb
Remarks/Alias	CD244; 2B4; NAIL; NKR2B4; Nmrk; SLAMF4; CD244 molecule
Species	Human
GeneID (Human)	51744
GeneID	51744
Immunogen	Recombinant protein corresponding to amino acids 22-221 of human CD244/2B4 (NP_001160135.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BZW8
Protein Weight	30kDa/36kDa/41kDa
Shipping	Ice bag