

Rabbit Anti-Human Cd244/2B4 Monoclonal Antibody For Flow Cytometry - Q9Bzw8

Item Number: CM0010537

Introduction

Rabbit monoclonal antibody targeting human CD244/2B4, tested for flow cytometry. Immunogen: recombinant protein (aa 22-221). Recognizes NK cell receptor, useful for innate/adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD244/2B4 mAb
Remarks/Alias	CD244; 2B4; NAIL; NKR2B4; Nmrk; SLAMF4; CD244 molecule
Species	Human
Gene ID (Human)	51744
Gene ID	51744
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-221 of human CD244/2B4 (NP_001160135.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BZW8
Protein Weight	30kDa/36kDa/41kDa
Shipping	Ice bag