

Rabbit Anti-Human Cd261/Trail-R1 Monoclonal Antibody For Flow Cytometry - 000220

Item Number: CM0010389



Introduction

-conjugated rabbit monoclonal antibody targeting human CD261/TRAIL-R1 (TNFRSF10A, DR4). Designed for flow cytometry (FC) applications. Reacts with human. Ideal for studying TRAIL-induced apoptosis and NF-kB activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD261/TRAIL-R1 mAb
Remarks/Alias	TNFRSF10A; APO2; CD261; DR4; TRAILR-1; TRAILR1; TNF receptor superfamily member 10a
Species	Human
Gene ID	8797 (Human)
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 124-232 of human CD261/TRAIL-R1 (NP_003835.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	000220
Protein Weight	50kDa
Shipping	Ice bag