

Rabbit Anti-Human Cd261/Trail-R1 Monoclonal Antibody For Flow Cytometry - O00220

Item Number: CM0010827



Introduction

High-quality -conjugated rabbit monoclonal antibody against human CD261/TRAIL-R1 (TNFRSF10A) for flow cytometry. Specifically designed for FC applications, this antibody targets the death receptor TRAIL-R1, enabling reliable detection of apoptotic pathways in human cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD261/TRAIL-R1 mAb
Remarks/Alias	TNFRSF10A; APO2; CD261; DR4; TRAILR-1; TRAILR1; TNF receptor superfamily member 10a
Species	Human
Gene ID (Human)	8797
Gene ID	8797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-232 of human CD261/TRAIL-R1 (NP_003835.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00220
Protein Weight	50kDa
Shipping	Ice bag