


Rabbit Anti-Human Cd262/Dr5/Trailr2 Monoclonal Antibody For Fc - O14763

Item Number: CM0006020



Introduction

The Rabbit anti-Human CD262/DR5/TRAILR2 Monoclonal Antibody is a fluorophore-conjugated antibody for flow cytometry (FC) detection of human TRAIL receptor 2 (DR5), a key mediator of apoptosis. From Rabbit, targets amino acids 56-182.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD262/DR5/TRAILR2 mAb
Remarks/Alias	TNFRSF10B; CD262; DR5; KILLER; KILLER/DR5; TRAIL-R2; TRAILR2; TRICK2; TRICK2A; TRICK2B; TRICKB; ZTNFR9; TNF receptor superfamily member 10b
Species	Human
GeneID (Human)	8795
GeneID	8795
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 56-182 of human CD262/DR5/TRAILR2 (NP_003833.4).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14763
Protein Weight	12kDa/45kDa/47kDa
Shipping	Ice bag