

Rabbit Anti-Human Cd263/Trail-R3 Monoclonal Antibody For Flow Cytometry - O14798

Item Number: CM0012133



Introduction

-conjugated rabbit monoclonal antibody against human CD263/TRAIL-R3 (TNFRSF10C), a decoy receptor inhibiting TRAIL-induced apoptosis. Validated for flow cytometry, ideal for apoptosis and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD263/TRAIL-R3 mAb
Remarks/Alias	TNFRSF10C; CD263; DCR1; DCR1-TNFR; LIT; TRAIL-R3; TRAILR3; TRID; TNF receptor superfamily member 10c
Species	human
GeneID (Human)	8794
GeneID	8794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-235 of human CD263/TRAIL-R3 (NP_003832.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14798
Protein Weight	27kDa
Shipping	Ice bag