

Rabbit Anti-Human Cd26/Dpp4 Mab For Flow Cytometry - P27487

Item Number: CM0004425



Introduction

Rabbit anti-Human CD26/DPP4 monoclonal antibody for flow cytometry. Recognizes human CD26, a T-cell activation antigen and dipeptidyl peptidase IV involved in immune regulation and cell adhesion. Recombinant rabbit mAb, UniProt P27487.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD26/DPP4 mAb
Remarks/Alias	CD26; ADABP; ADCP2; DPPIV; TP103
Species	Human
GeneID (human)	1803
GeneID	1803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-766 of human CD26/DPP4 (NP_001926.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag