

Rabbit Anti-Human Cd26/Dpp4 Monoclonal Antibody For Flow Cytometry - P27487

Item Number: CM0005478



Introduction

-conjugated rabbit monoclonal antibody targeting human CD26/DPP4. Validated for flow cytometry. Suitable for T-cell activation, dipeptidyl peptidase studies, and MERS-CoV receptor research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD26/DPP4 mAb
Remarks/Alias	DPP4; ADABP; ADCP2; CD26; DPPIV; TP103; dipeptidyl peptidase 4
Species	Human
GeneID (Human)	1803
GeneID	1803
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-766 of human CD26/DPP4(NP_001926.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P27487
Protein Weight	88kDa
Shipping	Ice bag