

Rabbit Anti-Human Cd27 Mab For Fc - P26842

Item Number: CM0010817



Introduction

-conjugated rabbit monoclonal antibody targeting human CD27 (TNFRSF7), validated for flow cytometry (FC). CD27 is a costimulatory immune-checkpoint receptor on T-cells, NK-cells, and B-cells, involved in T-cell activation and memory. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
GeneID (Human)	939
GeneID	939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-192 of human CD27 (NP_001233.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag