

Rabbit Anti-Human Cd27 Monoclonal Antibody For Fc - P26842

Item Number: CM0009753



Introduction

-conjugated rabbit monoclonal antibody targeting human CD27 (TNFRSF7) for flow cytometry. Immunogen: recombinant protein (aa 21-192). Suitable for T-cell, B-cell, and NK-cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
Gene ID (Human)	939
Gene ID	939
Immunogen	Recombinant protein (amino acids 21-192 of human CD27)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag