

Rabbit Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry - P26842

Item Number: CM0011684



Introduction

Rabbit monoclonal anti-human CD27 antibody (conjugate) validated for flow cytometry. Detects CD27/TNFRSF7, a costimulatory immune receptor on T, NK, and B cells. Cross-reacts with human CD27. Immunogen: recombinant human CD27 (aa 21-192).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
GeneID (Human)	939
GeneID	939
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-192 of human CD27 (NP_001233.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag