

Rabbit Anti-Human Cd270/Hvem Monoclonal Antibody For Flow Cytometry - Q92956

Item Number: CM0005704



Introduction

Rabbit monoclonal antibody against human CD270/HVEM (TNFRSF14), conjugated to . Validated for flow cytometry. Targets TNF receptor superfamily member 14 involved in immune cell signaling, herpesvirus entry, and immune checkpoint regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD270/HVEM mAb
Remarks/Alias	TNFRSF14; ATAR; CD270; HVEA; HVEM; LIGHTR; TR2; TNF receptor superfamily member 14
Species	Human
GeneID (Human)	8764
GeneID	8764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-162 of human CD270/HVEM (NP_003811.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q92956
Protein Weight	21kDa/30kDa
Shipping	Ice bag