


Rabbit Anti-Human Cd273/Pd-L2 Mab For Flow Cytometry - Q9Bq51

Item Number: CM0011814



Introduction

-conjugated rabbit monoclonal antibody targeting human PD-L2/CD273. Validated for flow cytometry. Suitable for studying immune checkpoint regulation and T-cell inhibition in cancer and autoimmune research.

[Learn More](#)

Parameter	Description
Product Name	Rabbit anti-Human CD273/PD-L2 mAb
Remarks/Alias	B7DC; Btdc; PDL2; CD273; PD-L2; PDCD1L2; bA574F11.2; PDCD1LG2
Species	Human
Gene ID (Human)	80380
Gene ID	80380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-219 of human CD273/PD-L2 (NP_079515.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BQ51
Protein Weight	21kDa/31kDa
Shipping	Ice bag