


Rabbit Anti-Human Cd28 Monoclonal Antibody For Fc - P10747

Item Number: CM0012036



Introduction

-conjugated rabbit monoclonal antibody against human CD28 (Tp44), designed for flow cytometry. This antibody targets the T-cell-specific surface glycoprotein CD28, enabling studies on T-cell activation, proliferation, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD28 mAb
Remarks/Alias	Tp44
Species	Human
GeneID (human)	940
GeneID	940
Immunogen	Recombinant protein of human CD28 Full length (P10747)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	P10747
Protein Weight	25kDa
Shipping	Ice bag