


Rabbit Anti-Human Cd28 Monoclonal Antibody For Flow Cytometry And T Cell Stimulation - P10747

Item Number: CM0010975



Introduction

Rabbit monoclonal antibody targeting human CD28 (Tp44), a T-cell costimulatory receptor. Suitable for flow cytometry and in vitro T-cell stimulation/activation. Ideal for immunology and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD28 monoclonal antibody
Remarks/Alias	Tp44
Species	Human
GeneID (Human)	940
GeneID	940
Immunogen	Recombinant protein of human CD28
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, in vitro T cell stimulation/activation
Cross Reactivity	Human
SWISS	P10747
Protein Weight	25kDa
Shipping	Ice bag