

Rabbit Anti-Human Cd298/Atp1B3 Mab For Fc - P54709

Item Number: CM0009067



Introduction

-conjugated rabbit monoclonal antibody against human CD298/ATP1B3, validated for flow cytometry. Targets the sodium/potassium-transporting ATPase beta-3 subunit, ideal for cell surface staining in IVD applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD298/ATP1B3 mAb
Remarks/Alias	CD298; ATPB-3
Species	Human
GeneID (human)	483
GeneID	483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-279 of human CD298/ATP1B3 (NP_001670.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P54709
Protein Weight	32kDa
Shipping	Ice bag