

Rabbit Anti-Human Cd298/Atp1B3 Monoclonal Antibody For Fc - P54709

Item Number: CM0008996



Introduction

-conjugated rabbit monoclonal antibody specifically recognizing human CD298/ATP1B3, the non-catalytic beta-3 subunit of the sodium/potassium-transporting ATPase. Suitable for flow cytometry (FC) with recombinant immunogen and validated cross-reactivity in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD298/ATP1B3 mAb
Remarks/Alias	CD298; ATPB-3
Species	Human
Gene ID (Human)	483
Gene ID	483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-279 of human CD298/ATP1B3 (NP_001670.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P54709
Protein Weight	32kDa
Shipping	Ice bag