


Rabbit Anti-Human Cd30 Monoclonal Antibody For Fc - P28908

Item Number: CM0010436



Introduction

-conjugated rabbit monoclonal antibody targeting human CD30 (Ki-1) for flow cytometry. Recognizes the 50-63kDa TNFR superfamily member involved in lymphocyte activation. Validated for human cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD30 mAb
Remarks/Alias	CD30; Ki-1; D1S166E
Species	Human
GeneID (Human)	943
GeneID	943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-379 of human CD30 (NP_001234.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28908
Protein Weight	14kDa/50kDa/63kDa
Shipping	Ice bag