


Rabbit Anti-Human Cd30 Monoclonal Antibody For Flow Cytometry - P28908

Item Number: CM0009893



Introduction

-conjugated rabbit monoclonal antibody targeting human CD30 (TNFRSF8, Ki-1). Recognizes extracellular domain (aa 19-379). For flow cytometry detection of CD30+ cells. Reactivity: human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD30 mAb
Remarks/Alias	CD30; Ki-1; D1S166E
Species	Human
Gene ID (Human)	943
Gene ID	943
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 19-379 of human CD30 (NP_001234.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28908
Protein Weight	14kDa/50kDa/63kDa
Shipping	Ice bag