


Rabbit Anti-Human Cd300A Monoclonal Antibody For Flow Cytometry - Q9Ugn4

Item Number: CM0008107



Introduction

-conjugated rabbit monoclonal antibody targeting human CD300a (CMRF35H, IRp60), validated for flow cytometry. Detects inhibitory receptor on NK cells, T cells, and mast cells. UniProt: Q9UGN4.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD300a mAb
Remarks/Alias	IRC1; IRC2; CLM-8; IRp60; IGSF12; CMRF35H; CMRF-35H; CMRF35-H; CMRF35H9; CMRF35-H9; IRC1/IRC2; CMRF-35-H9
Species	Human
Gene ID (Human)	11314
Gene ID	11314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-178 of human CD300a (NP_009192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UGN4
Protein Weight	33kDa/20kDa/12kDa/16.5kDa
Shipping	Ice bag