


Rabbit Anti-Human Cd300C Monoclonal Antibody For If/Icc, Fc - Q08708

Item Number: CM0011872



Introduction

-conjugated rabbit anti-human CD300c monoclonal antibody validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Detects CD300c, a cell membrane immune receptor expressed on monocytes, neutrophils, and lymphocyte subsets. For human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD300c mAb
Remarks/Alias	LIR; CLM-6; CMRF35; IGSF16; CMRF-35; CMRF35A; CMRF-35A; CMRF35A1; CMRF35-A1
Species	Human
GeneID (Human)	10871
GeneID	10871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-183 of human CD300C (NP_006669.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q08708
Protein Weight	25kDa
Shipping	Ice bag