


Rabbit Anti-Human Cd300F/Irem-1 Monoclonal Antibody For Flow Cytometry - Q8Tdq1

Item Number: CM0011478



Introduction

-conjugated rabbit monoclonal antibody targeting human CD300F/IREM-1 (CMRF35-like molecule 1), validated for flow cytometry. This inhibitory receptor regulates myeloid cell function, efferocytosis, and allergic responses. Ideal for immune cell profiling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD300F/IREM-1 mAb
Remarks/Alias	CD300LF; CD300f; CLM-1; CLM1; IREM-1; IREM1; IgSF13; LMIR3; NKIR; CMRF35-like molecule 1
Species	Human
GeneID (Human)	146722
GeneID	146722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-156 of human CD300F/IREM-1 (NP_620587.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ1
Protein Weight	18kDa/21kDa/26kDa/32kDa
Shipping	Ice bag