

Rabbit Anti-Human Cd301/Clec10A Monoclonal Antibody For Fc - Q8Iun9

Item Number: CM0010938



Introduction

-conjugated rabbit monoclonal antibody against human CD301/CLEC10A (CLEC10A), validated for flow cytometry. Recognizes a C-type lectin receptor on dendritic cells; specific for human. For immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD301/CLEC10A mAb
Remarks/Alias	HML, MGL, HML2, CD301, CLECSF13, CLECSF14, CD301/CLEC10A
Species	Human
GeneID (Human)	10462
GeneID	10462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-316 of human CD301/CLEC10A (NP_878910.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8IUN9
Protein Weight	28kDa/32kDa/35kDa
Shipping	Ice bag