

Abflo 488 Rabbit Anti-Human Cd301/Clec10A Monoclonal Antibody For Fc - Q8Iun9

Item Number: CM0004411



Introduction

ABflo 488-conjugated rabbit monoclonal antibody targeting human CD301/CLEC10A for flow cytometry. Detects the C-type lectin receptor involved in GalNAc recognition, immune modulation, and antigen presentation on myeloid cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD301/CLEC10A mAb
Remarks/Alias	HML; MGL; HML2; CD301; CLECSF13; CLECSF14; CD301/CLEC10A
Species	Human
Gene ID (Human)	10462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-316 of human CD301/CLEC10A (NP_878910.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8IUN9
Protein Weight	28kDa/32kDa/35kDa
Shipping	Ice bag