

Rabbit Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Fc - Q8Wtt0

Item Number: CM0007776



Introduction

-conjugated rabbit monoclonal antibody against human CD303/BDCA-2 (CLEC4C), validated for flow cytometry. Targets plasmacytoid dendritic cell marker. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD303/BDCA-2 mAb
Remarks/Alias	CLEC4C; BDCA2; CD303; CLECSF11; CLECSF7; DLEC; HECL; PRO34150; C-type lectin domain family 4 member C
Species	human
GeneID (human)	170482
GeneID	170482
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8WTT0
Protein Weight	21kDa/25kDa
Shipping	Ice bag