

Rabbit Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Flow Cytometry - Q8Wtt0

Item Number: CM0012287



Introduction

-conjugated rabbit monoclonal antibody against human CD303/BDCA-2 (CLEC4C) for flow cytometry. Recognizes human CD303. Suitable for plasmacytoid dendritic cell identification in research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD303/BDCA-2 mAb
Remarks/Alias	CLEC4C; BDCA2; CD303; CLECSF11; CLECSF7; DLEC; HECL; PRO34150; C-type lectin domain family 4 member C
Species	human
Gene ID	170482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-213 of human CD303/BDCA-2 (NP_569708.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8WTT0
Protein Weight	21kDa/25kDa
Shipping	Ice bag