

## Rabbit Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Flow Cytometry - Q8Wtt0

Item Number: CM0011860



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD303/BDCA-2 (CLEC4C). Validated for flow cytometry, ideal for immunophenotyping plasmacytoid dendritic cells. SWISS: Q8WTT0.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD303/BDCA-2 mAb                                                                                       |
| Remarks/Alias    | CLEC4C; BDCA2; CD303; CLECSF11; CLECSF7; DLEC; HECL; PRO34150; C-type lectin domain family 4 member C                    |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 170482                                                                                                                   |
| Gene ID          | 170482                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 45-213 of human CD303/BDCA-2 (NP_569708.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | FC                                                                                                                       |
| Cross Reactivity | Human                                                                                                                    |
| SWISS            | Q8WTT0                                                                                                                   |
| Protein Weight   | 21kDa/25kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                                  |