

Rabbit Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Flow Cytometry - Q8Wtt0

Item Number: CM0008854



Introduction

Monoclonal rabbit antibody against human CD303/BDCA-2, conjugated to , validated for flow cytometry. Ideal for studying plasmacytoid dendritic cells and antigen presentation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD303/BDCA-2 mAb
Remarks/Alias	DLEC; HECL; BDCA2; CD303; BDCA-2; CLECSF7; CLECSF11; PRO34150
Species	Human
GeneID (Human)	170482
GeneID	170482
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 45-213 of human CD303/BDCA-2 (NP_569708.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8WTT0
Protein Weight	21kDa/25kDa
Shipping	Ice bag