


Rabbit Anti-Human Cd304/Neuropilin-1 Monoclonal Antibody For Flow Cytometry - O14786

Item Number: CM0009792



Introduction

-conjugated rabbit monoclonal antibody targeting human CD304/Neuropilin-1 (SWISS O14786) for flow cytometry. Validated with recombinant immunogen encompassing the extracellular domain. Ideal for angiogenesis, neurodevelopment, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD304/Neuropilin-1 mAb
Remarks/Alias	NP1; NRP; BDCA4; CD304; VEGF165R
Species	Human
GeneID (Human)	8829
GeneID	8829
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-644 of human CD304/Neuropilin-1 (NP_003864.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14786
Protein Weight	68kDa/71kDa/103kDa
Shipping	Ice bag