

Rabbit Anti-Human Cd305/Lair1 Mab For Fc - Q6Gtx8

Item Number: CM0006117



Introduction

Rabbit anti-Human CD305/LAIR1 monoclonal antibody for flow cytometry (FC). Specific for human LAIR1, expressed on NK, T, B cells, monocytes, dendritic cells. Ideal for immune cell receptor analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD305/LAIR1 mAb
Remarks/Alias	LAIR1; CD305; LAIR-1; leukocyte-associated immunoglobulin-like receptor 1
Species	Human
Gene ID (human)	3903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-163 of human CD305/LAIR1 (NP_002278.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q6GTX8
Protein Weight	23kDa/29kDa/31kDa
Shipping	Ice bag