

Rabbit Anti-Human Cd307A/Fcrl1 Monoclonal Antibody For Fc - Q96La6

Item Number: CM0010267

Introduction

-conjugated rabbit monoclonal antibody targeting human CD307a/FCRL1 (Fc receptor-like protein 1). Validated for flow cytometry in human samples, supporting B-cell subset analysis and immunophenotyping research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD307a/FCRL1 mAb
Remarks/Alias	FCRH1; IFGP1; IRTA5; CD307a
Species	Human
Gene ID (Human)	115350
Gene ID	115350
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 17-303 of human CD307a/FCRL1 (NP_443170.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96LA6
Protein Weight	47kDa/40kDa/22kDa
Shipping	Ice bag