

Rabbit Anti-Human Cd307A/Fcrl1 Monoclonal Antibody For Fc - Q96La6

Item Number: CM0005122



Introduction

-conjugated rabbit monoclonal antibody targeting human CD307a/FCRL1 (FCRH1, IRTA5, CD307a). Validated for flow cytometry. Suitable for studying B-cell activation and signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD307a/FCRL1 mAb
Remarks/Alias	FCRH1; IFGP1; IRTA5; CD307a
Species	Human
Gene ID (Human)	115350
Gene ID	115350
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96LA6
Protein Weight	47kDa/40kDa/22kDa
Shipping	Ice bag