

Rabbit Anti-Human Cd307A/Fcrl1 Monoclonal Antibody For Flow Cytometry - Q96La6

Item Number: CM0004265

Introduction

-conjugated rabbit monoclonal antibody targeting human CD307a/FCRL1, validated for flow cytometry. Suitable for B-cell receptor signaling studies, immunophenotyping, and B-cell research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD307a/FCRL1 mAb
Remarks/Alias	FCRH1; IFGP1; IRTA5; CD307a
Species	Human
Gene ID (Human)	115350
Gene ID	115350
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 17-303 of human CD307a/FCRL1 (NP_443170.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96LA6
Protein Weight	47kDa/40kDa/22kDa
Shipping	Ice bag