

Rabbit Anti-Human Cd307B (Fcrl2) Monoclonal Antibody For Fc - Q96La5

Item Number: CM0004347

Introduction

-conjugated rabbit monoclonal antibody targeting human CD307b (FCRL2) for flow cytometry. Validated on human samples. CD307b may regulate B cell development and is preferentially expressed in memory B cells.

[Learn More](#)

Product Name	Rabbit anti-Human CD307b mAb
Remarks/Alias	FCRH2, IFGP4, IRTA4, SPAP1, CD307b, SPAP1A, SPAP1B, SPAP1C
Species	Human
Gene ID	79368
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 20-395 of human CD307b (NP_110391.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	Q96LA5
Protein Weight	56kDa
Shipping	Ice bag