

Rabbit Anti-Human Cd307B Monoclonal Antibody For Fc - Q96La5

Item Number: CM0005869

Introduction

-conjugated rabbit monoclonal antibody against human CD307b (FCRL2/IRTA4). Validated for flow cytometry, suitable for B cell development and memory B cell research.

[Learn More](#)

Specification	Details
Product Name	Rabbit anti-Human CD307b mAb
Remarks/Alias	FCRH2; IFGP4; IRTA4; SPAP1; CD307b; SPAP1A; SPAP1B; SPAP1C
Species	Human
GeneID (Human)	79368
GeneID	79368
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-395 of human CD307b (NP_110391.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96LA5
Protein Weight	56kDa
Shipping	Ice bag