

Rabbit Anti-Human Cd307B Monoclonal Antibody For Flow Cytometry - Q96La5

Item Number: CM0005863

Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD307b (FCRL2) for flow cytometry. Validated for human cross-reactivity, 56kDa protein weight, recombinant immunogen. Ideal for B-cell immunology research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD307b mAb
Remarks/Alias	FCRH2; IFGP4; IRTA4; SPAP1; CD307b; SPAP1A; SPAP1B; SPAP1C
Species	Human
GeneID	79368
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-395 of human CD307b (NP_110391.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96LA5
Protein Weight	56kDa
Shipping	Ice bag