


Rabbit Anti-Human Cd314 Monoclonal Antibody For Fc - P26718

Item Number: CM0010636



Introduction

-conjugated rabbit monoclonal antibody against human CD314 (NKG2D), validated for flow cytometry. Enables direct detection of activating NK cell receptor expression in immunological and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD314 mAb
Remarks/Alias	KLR; CD314; NKG2D; NKG2-D; D12S2489E
Species	Human
Gene ID (Human)	22914
Gene ID	22914
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26718
Protein Weight	25kDa
Shipping	Ice bag