

Rabbit Anti-Human Cd314 (Nkg2D) Monoclonal Antibody For Fc - P26718

Item Number: CM0008716



Introduction

-conjugated rabbit monoclonal antibody against human CD314 (NKG2D), validated for flow cytometry. Ideal for immunosurveillance studies and NK/T-cell activation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD314 mAb
Remarks/Alias	KLR; CD314; NKG2D; NKG2-D; D12S2489E
Species	Human
Gene ID (Human)	22914
Gene ID	22914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-216 of human CD314 (NP_031386.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26718
Protein Weight	25kDa
Shipping	Ice bag