

Rabbit Anti-Human Cd317/Bst2 Mab For Fc - Q10589

Item Number: CM0009447



Introduction

-conjugated rabbit anti-human CD317/BST2 (Tetherin/HM1.24) monoclonal antibody, validated for flow cytometry (FC). Targets a 20 kDa antiviral restriction factor; suitable for innate immunity and viral budding research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD317/BST2 mAb
Remarks/Alias	CD317; HM1.24; TETHERIN
Species	Human
Gene ID (Human)	684
Gene ID	684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-160 of human CD317/BST2 (NP_004326.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q10589
Protein Weight	20kDa
Shipping	Ice bag