

Rabbit Anti-Human Cd317/Bst2 Monoclonal Antibody For Fc - Q10589

Item Number: CM0012269



Introduction

-conjugated rabbit monoclonal antibody targeting human CD317/BST2 (Tetherin), validated for flow cytometry (FC). Ideal for studying antiviral restriction and immune signaling in human samples. 18-19 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD317/BST2 mAb
Remarks/Alias	CD317; HM1.24; TETHERIN
Species	Human
Gene ID (Human)	684
Gene ID	684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-160 of human CD317/BST2 (NP_004326.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q10589
Protein Weight	18kDa/19kDa
Shipping	Ice bag